



Sampling for Laboratory Testing

	Skin lesions	Feather follicles	Cloacal swabs	Brain	Trachea/swabs	Lungs	Liver	Spleen	Pancreas	Kidneys	Oviduct	Thymus	Bone marrow	Bursa fabricius	Proventriculus/ Gizzard	Caecal tonsils	various organs (tumors) Nerves	Joints	Virus culture possible *
AE			✓✓					✓											
AI		✓✓			✓✓	✓✓		✓							✓✓				X
AmPV (ART)					✓✓	✓✓									✓				
CAV								✓✓				✓✓							
EDS		✓									✓✓								
FAdV		✓					✓✓	✓	✓✓					✓	✓✓				X
Gumboro (IBD)		✓											✓✓	✓					X
IBV		✓✓			✓✓	✓✓				✓	✓				✓✓				X
ILT					✓✓	✓									✓				X
Marek's/ Leukosis	✓✓							✓✓						✓✓		✓✓	✓✓		
Mycoplasma		✓✓			✓✓						✓✓				✓			✓✓	
ND		✓✓			✓✓	✓✓									✓✓				X
Pox	✓✓				✓✓														X
Reo		✓					✓	✓	✓✓						✓✓			✓✓	X

✓✓ most preferred material

✓ also possible

*Isolation of pathogens require specific transport medium and transport conditions