



# Probenmaterial für Labortests

	Hautverletzungen	Federfölikel	Kloakentupfer	Gehirn	Trachealtupfer	Lunge	Leber	Milz	Bauchspeicheldrüse	Nieren	Eileiter	Thymus	Knochenmark	Schleimbeutel	Drüsen-/Muskelmagen	Ceacal-Tonsillen	versch. Organe (Tumore)	Nerven	Gelenke	Virusanzucht möglich*
AE				✓✓					✓											
AI			✓✓		✓✓	✓✓		✓							✓✓					X
AmpV (ART)					✓✓	✓✓									✓					
CAV								✓✓			✓✓	✓✓								
EDS			✓							✓✓					✓					
FAdV			✓				✓✓	✓	✓✓						✓	✓✓				X
Gumboro (IBD)			✓										✓✓	✓						X
IBV			✓✓		✓✓	✓✓				✓	✓				✓✓					X
ILT					✓✓	✓									✓					X
Marek's/Leukosis		✓✓						✓✓							✓✓		✓✓	✓✓		
Mycoplasma			✓✓		✓✓					✓✓					✓				✓✓	
ND			✓✓		✓✓	✓✓									✓✓					X
Pox	✓✓				✓✓															X
Reo			✓				✓	✓	✓✓						✓✓			✓✓		X
✓✓	bevorzugtes Material																			
✓	alternatives Material																			

\*Die Isolierung von Krankheitserregern erfordert ein spezifisches Transportmedium und spezifische Transportbedingungen